

# URIC ACID (DST)

(URICASE METHOD)

ENZOPAK

Last update 02-2026

**Ref.** CC2-URK.20N, 5x25 ml  
CC2-URK.20NU, 5x60 ml

## INTENDED USE

Reagent kit for quantitative estimation of Uric acid in serum or plasma.

## PRODUCT HIGHLIGHTS

- Reagent free from interference due to reducing substances
- Specially related chromogen for increased sensitivity.

## INTRODUCTION

Uric acid is produced by the action of xanthine oxidase on xanthine and hypoxanthine, which are products of nucleic acid degradation. The chemical method of estimation of uric acid is based on the uric acid forming a blue phosphotungstate complex. However, this method in serum or plasma continue to produce colour by phosphotungstate complex. In 1980, Fossati et, al described a modified Trinder method for the assay of uric acid.

ENZOPAK Uric Acid is formulated on this method with the advantage that the method is specific and sensitive.

## DIAGNOSTIC SIGNIFICANCE


Uric Acid levels in any subject may vary during a day from time to time and from one to another day. The changes represent the contents in diet and metabolism besides pathological condition.

The diagnostic values showing increased uric acid levels are in kidneys failure and gout. Elevated levels are found in acute infectious diseases, severe uremia, toxemia of pregnancy, leukemia and some malignant diseases or conditions.

The decreased levels are found due to drugs and hormone treatment like adrenocorticotrophic hormone (ACTH).

## PRINCIPLE

Uricase is a very specific enzyme acting on Uric Acid, end products being allantoin and peroxide. Peroxidase is used to utilize hydrogen peroxide (proportional to uric acid concentration) to convert chromogen to coloured complex. The intensity of colour produced is proportional to uric acid concentration and is measured photometrically at 520 nm (500-550 nm or GREEN filter).



TBHB : 2,4,6-Tribromo-3-hydroxy benzoic acid.

## PRESENTATION

All reagents to be stored at 2-8°C

	No. of Bottles	
	5x25 ml	5x60 ml
• 1 Uric Acid (Enzyme/Chromogen)	5	5
• 2 Uric Acid (Buffer)	1	5
• Uric Acid Standard (6 mg/dl)	1	1

## FINAL REAGENT COMPOSITION

Active Ingredients	Concentration
• Uricase	≥ 200 U/L
• POD	≥ 1500 U/L
• 4-AAP	0.2 mmol/L
• TBHB	0.5 mmol/L
• Buffer	100 mmol/L
• Detergent	1 mmol/L

pH 8.2 ± 0.5 at 25°C

Uric Acid Standard (6mg/dl)

Also contains non-reactive fillers and Stabilizers.

## PRECAUTION

ENZOPAK Uric acid is for *IN-VITRO* diagnostic use only.

Reagent contains Sodium Azide. DO NOT INGEST.

## PREPARATION OF WORKING REAGENT

### For 5 x 25 ml

Dissolve the contents of one vial of 1 Uric Acid (Powder) with 25 ml of 2 Uric Acid (Buffer). Mix well to dissolve, wait for 5-minutes before use.

### For 5 x 60 ml

Carefully transfer the content of 1 Uric Acid (Powder) into the bottle containing 60ml of 2 Uric Acid Buffer. Mix well to dissolve. Wait for 5 minutes before use.

## REAGENT STORAGE & STABILITY

All reagents included in the kit are stable at 2-8°C until the expiry date indicated on the label. Working reagent is stable for 6 months at 2-8°C. A slight rise in the blank reading due to stored life will not affect the test results.

## SPECIMEN COLLECTION

Fresh, clear fasting serum with no hemolysis is the specimen of choice. However, plasma collected from blood using heparin as an anticoagulant may be used.

## REACTION PARAMETERS

• Type of Reaction	:	End Point
• Wavelength	:	520 nm (500-550 nm)
• Flowcell Temperature	:	37°C
• Incubation	:	5 min. at 37°C
• Std. Concentration	:	6 mg/dl
• Sample Volume	:	20 µl (0.02 ml)
• Reagent Volume	:	1.0 ml
• Light Path	:	1.0 cm
• Zero setting with	:	Reagent Blank

## TEST PROCEDURE

Pipette Into Test Tubes	Blank	Standard	Test
Working Reagent (ml)	1.0	1.0	1.0
Standard (ml)	-	0.02	-
Sample (ml)	-	-	0.02

Mix and incubate at 37°C for 5 minutes and read absorbance of test and standard against reagent blank at 520 nm (500-550 nm or GREEN filter).

## STABILITY OF FINAL COLOUR DEVELOPED

The colour of reaction mixture is stable for one hour when protected from direct light.

# URIC ACID (DST)

(URICASE METHOD)

**ENZOPAK**

Last update 02-2026

## TEST RESULTS

$$\text{Uric Acid conc. (mg/dl)} = \frac{\text{Absorbance of Test}}{\text{Absorbance of Std.}} \times 6$$

To convert Uric Acid concentration (mg/dl) to micromoles/liter use the following equation.

$$\text{Micromoles/lit.} = \text{mg/dl} \times 59.5$$

## NORMAL VALUES

Male : 3.2-7.0 mg/dl (190.4 – 416.5  $\mu\text{mol/l}$ )

Female : 2.6-6.0 mg/dl (154.7 – 357  $\mu\text{mol/l}$ )

## LINEARITY

This method is linear upto 20 mg/dl. For sample values higher than linearity limit, dilute the samples suitably with 0.9 % saline and repeat the assay. Apply proper dilution factor to calculate the test results by applying dilution factors.

## REFERENCES

TRIVEDI, R.C., REBBAR, L.BERKA, E., AND Strong, I: Clin. Chem 24:1908 (1978).

FOSSATI, P., PRINCIPLE, L., BERTI, G: Clin Chem 26, 227-231 (1980).



Regd. Office.3/7, B.I.D.C., Gorwa, Vadodara 390 016 (INDIA)

Web: [www.reckondiagnosics.com](http://www.reckondiagnosics.com) Ph: +91-265-2281631

Email: [mail@reckondiagnosics.com](mailto:mail@reckondiagnosics.com)

An ISO 13485:2016 & GMP Certified Company

